

MIX 3

The **MIX-3** is a manual device for mixing two or more products which when shaken produce an emulsion. The MIX-3 is suitable, for example, for mixing gel with isopentane or other types of low-boiling-point liquids or liquefied gases; or for different liquids which when shaken produce a single final product.

The mixer is suitable for use in laboratory testing procedures and/or small production lines.

The unit comprises a base which holds a funnel, a two-piston cylinder with a pyrex® chamber, and the controls. The MIX-3 is very easy to use: the products to be mixed are poured through the funnel into the pyrex® mixing chamber.

A switch activates the cylinder that mixes the products. Once the desired degree of mixing has been achieved, the mix can be injected into a can using an injector.

